SCIENCE
CAN

MICROSCOPE VISION MANUAL

SCIENCE
CAN

VN:1602312025ig0001

TOPBRIGHT ANIMATION CORPORATION

TEL: +86 0571 86879395 Address: 6/F, The Agriculture building, Anyang district, Ruian, Zhejiang, China

TOPBRIGHT GmbH Brunnengasse 65, 90402 Nürnberg, Germany

E-mail: contact@topbrighttoys.com Website: <http://www.topbrighttoys.com>

Manufacturer: Zhejiang Topbright Toys Co., Ltd.

Address: No.2, Yangfan Road, Bailongshan Street, Yunhe County, Lishui City, Zhejiang Province, China

Sea&Mew Accounting Ltd Electric Avenue Vision 25, London, Enfield EN3 7GD

©2025 TOPBRIGHT ANIMATION all rights reserved. MADE IN CHINA

8-14
years

STEAM
EDUCATIONAL PRODUCT

WARNING CHOKING HAZARD - SMALL
PARTS, NOT SUITABLE FOR CHILDREN UNDER 3 YEARS.

CONTENT/INHALT/CONTENU /CONTENIDO/INHOUD

EN

THE MICROSCOPE VISION	P. 01
THE BUTTONS	P. 02
BASIC KNOWLEDGE	P. 03
HOW TO USE THE MICROSCOPE	P. 06

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Dive into the exciting world of the microcosm!

Have you ever wanted to know what a plant fiber looks like? Or see what is hidden in water? With a microscope, we see things that our eyes can't otherwise see. It magnifies these things and helps us look at them greatly magnified.

The history of the microscope

The history of the microscope begins a long time ago. In the 17th century, a clever man named Anton van Leeuwenhoek built one of the first microscopes. At the time, no one knew there were so many tiny things to discover. With his microscope, Leeuwenhoek saw tiny animals in water samples, opening up a new world.

But how does a microscope work?

A microscope has two parts: the eyepiece (through which you look) and the objective (which magnifies things). When you look through the eyepiece, the lens takes tiny things and makes them bigger so you can see them. It's like you're immersed in a miniature world!

And what do we use microscopes for?

Microscopes help scientists to discover tiny details. They are used in medicine to study diseases, in biology to study cells, and in chemistry to study tiny particles. But you can also have fun with them by looking closely at leaves, hair, or even insects.


Microscope with screen

(The microscope can be used independently.)

Lab base with object stage

Round specimen slide with sample x5

Round empty specimen slide x2

Rectangular empty specimen slide x2

Test tube x1

Tweezers x1

Sample vial x1

Pipettex1

USB C Cable x1

Petri dish x1

SD card x1

The Button	Function	Picture
Power Button	Power On/Off: Press and hold for 1 second.	
Shoot/Confirm Button	Take Photo: In the shooting interface, press briefly to take a photo.	
	Record Video: In the shooting interface, press and hold for 1 second to start recording, press briefly again to stop recording.	
	Video Playback: In the gallery interface, press briefly to play the video, press briefly again to pause.	

THE BUTTONS

The Button	Function	Picture
Gallery Button	Enter Gallery: In the shooting interface, press briefly to enter the gallery interface.	
	Delete Photo/Video: In the gallery interface, press and hold for 1 second to display the delete interface, then use the confirm button to choose whether to delete.	
	Return to Shooting Mode: In the gallery interface, press briefly to return to the shooting interface.	
Direction/ Magnification Button (Left)	Decrease Magnification: In the shooting interface, press briefly to decrease the magnification (1x is the default magnification).	

THE BUTTONS

The Button	Function	Picture
Direction/ Magnification Button (Left)	Switch Gallery: In the gallery interface, press briefly to switch to the previous photo/video.	
Direction/ Magnification Button (Right)	Increase Magnification: In the shooting interface, press briefly to increase the magnification (4x is the maximum magnification).	
	Switch Gallery: In the gallery interface, press briefly to switch to the next photo/video.	

Question 1: What do I do if the image is not in the center of the field of view?

1. Move your object of observation.
2. If the specimen is not centered (front, back, left or right) on the screen (as shown in Figure 1), move the specimen in the appropriate direction (front, back, left or right) until the image on the screen is aligned with the position of the specimen (as shown in Figure 2).

Figure 1

Figure 2

Question 2: What do I do if the image in the field of view is blurred?

Handheld microscope without object stage:

1. Make sure that the bottom of the microscope is close to the slide (as shown in Figure 3).
2. Turn the focusing wheel and adjust the image until it is in focus (as shown in Figure 4).
3. Lightly press the magnification button to select the desired magnification (as shown in Figure 5).

Figure 3

Figure 4

Figure 5

Microscope attached to the object stage:

1. Make sure that the microscope is in the lowest position.
2. Turn the focusing wheel and adjust the image until it is in focus (as shown in Figure 6).
3. Lightly press the magnification button to select the desired magnification (as shown in Figure 7).


Figure 6

Figure 7

How to save and transfer your photos and videos to a computer or smartphone

*Data transfer and computer display are supported on some Windows, macOS and Android systems (iOS and iPadOS are not supported).

1. Switch off the microscope first. Then connect the microscope to a computer or Android phone/tablet that supports OTG via the USB-C data cable.
2. Switch on the microscope and then briefly press the direction button to select the menu options
3. Briefly press the confirmation button to confirm.
4. Open the "Camera" or other application that can access the computer's camera.
5. Click on the "Flip camera" button
3. Press the photo/video button to record your observations.

The above steps refer to the Windows system. For other systems, you may need to download an application that supports USB or OTG cameras and follow the device-specific instructions.

Note: Some Windows, macOS and Android systems may not be able to view photos and videos stored on the user-supplied memory card. For security reasons, devices other than this product cannot access the photos stored on the internal memory.

1. Portable handheld mode

Required accessories: Microscope, specimen slide


Step 1: Press and hold the microscope power button for 1 second to turn on the microscope. The screen will display the imaging surface. Aim the objective lens at the area of the object you want to observe.

Step 2: If you are using a specimen slide, move the slide back and forth until the contents of the observation object appear on the microscope screen. (Note: The target should be directly under the center of the objective lens).

Step 3: Adjust the focusing wheel at the bottom of the microscope by hand until the image is in focus, then press the "Magnification" button on the microscope screen to select the desired magnification.

Note: The microscope slide can be replaced with other observation objects (such as leaves, etc.), but the steps remain the same.

2. The Laboratory Station

3. Laboratory Bench Mode

Required accessories: Object stage, experimental tools (tweezers, pipette, sample vessel, test tube)

Step 1: Place the tools in the left and right tool trays of the lab bench.

Step 2: Pull open the left drawer of the lab bench and take out the Petri dish and the round and rectangular specimen slides to prepare the experiment.

4. Rotatable object stage

Required accessories: Microscope, object stage, microscope slide

Step 1:

Assemble the microscope and stage as shown in the illustration. Turn the stand lifting wheel to raise the microscope to the highest position and then switch on the main unit.

Step 2:

Round specimen slide: Place the round specimen slide in the round area on the stage and secure it.

Rectangular specimen slide: Place the rectangular specimen slide in the rectangular specimen area on the stage and secure it.


Round specimen slide


Rectangular specimen slide

4. Rotatable object stage

Step 3: Rotate the stage to place the slide under the microscope.
Turn and adjust the “stand lifting wheel” by hand to bring the microscope to the lowest position.

Step 4:
Adjust the focus wheel by hand until the image is in focus.
Then press the “Magnification” button on the microscope screen to select the desired magnification.

What your laboratory equipment can do:

Required accessories: tools, observation objects (such as leaves)

1. Place the specimens that you have prepared and want to look at under the microscope on the empty slides. Seal each specimen between the two parts of the slide to keep them in the condition you need for your observations.

2. How to use the rectangular blank specimen slide:

Lay the slide down, put your specimen in the circle in the center and close the slide with the lid as shown in the illustration (the specimen should be small and thin enough for examination in the slide).

3. Use the test tube to measure the liquids for your experiments.

4. Use the pipette if you want to apply drops of water to the slide for your experiments.

What your laboratory equipment can do:

5. Use the petri dish to store your observation objects.

6. Use the tweezers and the sample vessel to search for new observation objects.

7. The tweezers will also help you to place your specimens on a specimen slide.

Note: Make sure that the microscope housing does not come into contact with water during use or cleaning.

WARNING!

Not suitable for children under 3 years.

Small parts, choking hazard.

Please check the list of equipment in advance to ensure that the correct parts are included in the box. Keep the packaging and instructions as they contain important information!

ATTENTION!

Any packaging material, including adhesive tapes, plastic pieces, binding wires, tags, etc, is not a part of the product and is not safe for children to play with. For the safety of your children, please dispose the above packaging materials properly. Please keep the contact details for further enquiries.

Connecting the charging cable:

1. Locate the connection for the USB C cable on the back of the main device and plug in the USB C charging cable.
2. Connect the charging cable to the power socket.
3. Charge the microscope when the battery is low. The estimated charging time is 3 hours.

Pre-Use Checks:

1. Check if the body of the product is damaged.
2. Ensure the USB charging cable is not damaged.

For children, use under adult supervision.

This product does not come with a power adapter. If using an adapter for power supply, make sure it meets CCC certification standards.

Maintenance:

1. Wipe the toy housing gently with a piece of soft, dry cloth to keep the frame clean.
2. Keep the main unit away from sunlight and any other heat source.
3. Please remove the batteries if the toy is to be left unused for a long time.
4. Do not touch the main unit with any hard object or try to disassemble it.
5. Do not put the main unit in water; please keep it from moisture.

Cautionary notes on the use of lithium batteries:

1. There is a risk of fire, explosion, and combustion. Do not disassemble, crush, heat above 60°C (140°F), or burn. Avoid using near heat sources like fire or heaters.
2. Do not submerge the product in water. Doing so could damage the internal protection circuit and function, or cause abnormal chemical reactions that may lead to overheating, smoking, deformation, or combustion.
3. This product and battery must be recycled by a professional recycling center. Do not dispose of them in household garbage.
4. The charging cable is not a toy and should be kept away from children. It should only be used with the provided device and should not be connected to more than the recommended number of power sources.
5. Regularly check the charging cable for damage to the power line, plug, outer casing, or other components. If damage is found, the charging cable should be discarded.
6. Always use the recommended adapter for charging the product. Overloaded adapters can cause the battery to overheat or become damaged. Do not charge the product for more than 24 hours.
7. Disconnect the toy from the adapter before cleaning.
8. The toy should only be connected to the specified devices.
9. This toy is not suitable for children under three years old.

Common Troubles and Troubleshooting:

If the toy stops responding for some reason, please:

1. Shut down the main unit, wait for a few minutes, and then turn it back on before using it again.
2. If the power cannot be turned off, locate the reset hole on the back of the main unit and briefly press it to reset. After restarting, the main unit will reboot.
3. Locate the reset hole on the back of the main unit, briefly press it to reset, and upon restarting, the main unit will restore to factory settings.

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